



Annotation Tool User Guide

- Click on **annotation.html** which is present on the desktop.
- After clicking, it should open in chrome




- Click on the Choose Files.
 - File explorer should open up, navigate to the Dicom directory where volume is present, select all the images within this folder. You can use a (Ctrl+a) shortcut to select all the files.



Click on Open


After clicking on Open, all the Dicom images selected all shown below.



- After annotating images with the help of radio buttons, click on 'Retrieve Annotation as CSV'. A popup window will appear.

- We need to provide a folder location where Dicom volume was present, this is required by the user because javascript doesn't allow the capture of file location when we are uploading the files. Also, provide an annotation file name.

- Click on the Download button, after clicking a CSV file will be generated by the tool which will contain the mapping of the Dicom 2d file and category it is assigned to by the user. The file can be found in the downloaded folder.



For better organization of these annotation files, please save this file the annotation file to the corresponding volume directory. The tool only creates a file and save it to the download folder.